Standard Operating Procedures # 12

Lab Name: FM Global Fire Phenomenon L	_aboratory		
Lab Location: Building: <u>J.M. Patterson (083</u>) Room		Room #:	3225
Laboratory director: <u>Michael J. Gollner</u>			Page _1_of _1_
1 Process	Operation of p	ilsed Nd-YA	Glaser (532 nm)

1. Process:	Operation of pulsed Nd-YAG laser (532 nm) and water tank.	
2. Hazards:	Micro-particle/fluorescent dye: PSL, Rh-6G. Facility: pulsed laser.	
3. Personal Protective Equipment:	Safety goggles and gloves with dye. Safety goggles, lattice gloves and dust masks with micro-particles. Laser safety goggles while laser is in operation.	
 Engineering \ Ventilation Controls: 	Not applicable.	
5. Special Handling Procedures Storage Requirements:	Laser safety barrier/curtains must be closed all around the water tank while the laser is in use.	
 Spill containment/ Accident Procedures: 	Stop work if any Rh-6G spills, wash contaminated parts with fresh water and notify the lab manager.	
7. Waste Disposal	All residual water+Rh-6G or water+PSL micro-particles should be placed in a special metal can and disposed of according to environmental safety procedures.	
8. Required Approvals:	Lab workers must be trained by lab manager in proper use of the Rh-6G dye and PSL micro-particles.	
	Campus laser safety training should be performed before working with the laser	
9. Decontamination:	No decontamination needed.	
10. Designated Areas:	Water tank.	